

RingJect™

Preloaded Disposable CTR Injector

97.8% IOL centration
at 12 months postoperative¹




RingJect™ Preloaded Capsular Tension Ring (CTR)

- Easy to use, single-use injector
- Implant clockwise or counterclockwise
- Small incision: 2.2 mm
- Ready to use in a sterile blister

Capsular Tension Rings are indicated for the stabilization of weakened, broken or missing zonules¹.

RingJect™

Preloaded Disposable CTR Injector




CTR is made of unique compression-molded PMMA: Extremely flexible



The ends of the rings are gradually formed to “tip up” like a ski tip – this allows the CTR to be easily guided in the capsular bag


Operating Instructions



1 Remove the RingJect™ from its blister pack



2 Depress the plunger to engage the CTR; stop pushing after the “click”



3 Slowly release the plunger; the CTR will load in the injector

4 The injector is now loaded and ready for insertion into the eye



IMPLANTATION



Counterclockwise implantation



Clockwise implantation

One instrument for clockwise or counterclockwise implantation. The arrow on the base of the tip of the injector shows the direction in which the capsular tension ring will exit the injector.

Features & Benefits

Preloaded: Self-loading, single packaging is a valuable addition to your surgical armamentarium

Injector: Easy to use

Total control with an action/retraction mechanism

Implantation – clockwise or counterclockwise indicators on the injector

Long, small tip is especially convenient in cases of deep-set eyes

Beveled tip for easy entry into the incision

Capsular Tension Ring: Circular expansion and stabilization of the capsular bag

Safe IOL centration in eyes with zonular dehiscence

Stabilized conditions during phacoemulsification surgery

Compression-molded PMMA: extremely flexible and strong

“Ski-tip” feature at the end of the ring allows the CTR to be easily guided in the capsular bag

RingJect™ Model 375 BT 1G00

Compression: From 12 mm to 10 mm
Material: PMMA | Flexible Ring
Available in: Clear PMMA
Overall Ø: 12 mm

RingJect™ Model 376 BT 1G00

Compression: From 13 mm to 11 mm
Material: PMMA | Flexible Ring
Available in: Clear PMMA
Overall Ø: 13 mm