


# TECNIS Simplicity™ — The Next-Generation TECNIS® IOL Disposable Delivery System


Simplifying delivery of the IOL you prefer.



## BENEFITS:




Easy 3-step process -  
hydrate, advance, deliver



Advanced design for a  
smooth and controlled  
IOL delivery<sup>2</sup>



Flexibility to use BSS  
or OVD for hydration<sup>2</sup>



For more information on TECNIS Simplicity™ Delivery System,  
please contact your Johnson & Johnson Vision Representative.

**TECNIS** SIMPLICITY™  
DELIVERY SYSTEM

**Johnson & Johnson** VISION

# The Next-Generation TECNIS® IOL Disposable Delivery System



**TECNIS  
EYHANCE™**

<b>Product Code:</b>	<b>DIB00</b>	
<b>OPTIC CHARACTERISTICS</b>		
<b>Power Range:</b>	+5.0 D to +34.0 D in 0.5 diopter increments	
<b>Diameter:</b>	6.0 mm	
<b>Shape:</b>	Biconvex, anterior aspheric surface	
<b>Material:</b>	UV-blocking hydrophobic acrylic	
<b>Refractive Index:</b>	1.47 at 35°C	
<b>Edge Design:</b>	ProTEC 360° Edge; frosted edge design; continuous uninterrupted posterior square edge	
<b>OPTICAL BIOMETRY/APPLANATION ULTRASOUND BIOMETRY</b>		
<b>A-constant:</b>	119.3*	118.8 <sup>†</sup>
<b>Theoretical AC Depth<sup>1</sup>:</b>	5.7 mm*	5.4 mm
<b>Surgeon Factor<sup>1</sup>:</b>	1.96 mm*	1.68 mm
<b>HAPTIC CHARACTERISTICS</b>		
<b>Lens Overall Diameter:</b>	13.0 mm	
<b>Style:</b>	C	
<b>Design:</b>	Tri-Fix, haptics offset from optics, 1-piece lens	
<b>Incision size:</b>	2.2 mm - 2.4 mm	

\*Derived from clinical evaluation results of the Tecnis® 1-Piece IOL Platform for optical biometry.

<sup>†</sup>A-Constant theoretically derived for ultrasound biometry.

**References:** **1.** Holladay, J. T. International Intraocular Lens & Implant Registry 2004. J Cataract Refract Surg 30, 207-29 (2004). **2.** TECNIS EYHANCE IOL with SmartLOAD™ Delivery Technology [DIB00 Package Insert] Z311246P Rev. B Revision date 08/2019. Santa Ana, Calif. Johnson & Johnson Surgical Vision, Inc.

**For Healthcare Professionals Only. Please read the Directions for Use and consult our specialists if you have any questions.**