




J&J Vision: Commonly Used Phaco Tips and Sleeves Guide

Phaco tips	INCISION SIZE: 2.2 mm to 2.4 mm	INCISION SIZE: 2.4 mm to 2.7 mm	INCISION SIZE: 2.7 mm to 3.0 mm
21 Gauge	Light Blue Sleeve (OPOS21LD) <i>Single Use</i> 	-	-
20 Gauge	-	Yellow Sleeve * (OPOS20L) <i>Reusable</i> 	-
19 Gauge	-	-	Blue Sleeve (OPOS19L) <i>Reusable</i> 

* For 2.4 mm incision size only




De Ceunynck
OPHTHALMOLOGY


Johnson & Johnson VISION




0° Bevel



15° Bevel



30° Bevel



45° Bevel



Straight Tip



Bent Tip
20° Angle

Most commonly used Phaco tips in EMEA

Straight tips	19 Gauge	20 Gauge	21 Gauge
0° Bevel	OPOR0019L, <i>reusable</i>		
15° Bevel	OPOR1519L, <i>reusable</i>	OPOR1520L, <i>reusable</i>	OPOR1521G, <i>disposable</i>
30° Bevel	OPOD3019L, <i>disposable</i>	OPOD3020L, <i>disposable</i>	OPOR3021G, <i>disposable</i>
	OPOF3019L, <i>disposable</i>	OPOR3020L, <i>reusable</i>	OPOR3021R, <i>reusable</i>
	OPOR3019L, <i>reusable</i>		
Bent tips			
15° Bevel			
30° Bevel		OPOCR3020L, <i>disposable</i>	OPOCR3021L, <i>disposable</i>
45° Bevel			

