

HEALON EndoCoat™ PRO

Enhances eye safety in cataract procedures

Offering premium protection for patient's eyes, but now with:^{1, 2}

- Animal-free Hyaluronic Acid¹
- An easy-to-mount cannula guard with a side-mounted snap-on design for enhanced ease of use^{1, 2}
- An ergonomic delivery system designed for ease of use throughout surgery^{1, 2}



Safeguarding ocular integrity during procedures with HEALON EndoCoat™ PRO dispersive OVD*:

- Protects the corneal endothelium and other ocular tissue^{1, 2}
- Enhances visualization during the surgical procedure²
- Maintains a deep chamber during anterior segment surgery^{1, 2}

*Ocular integrity refers to protection of corneal endothelium

 **Healon™**
FAMILY of OVDs

FROM THE MAKERS OF **TECNIS™**

Johnson & Johnson

HEALON EndoCoat™ PRO

3% Sodium Hyaluronate

SPEC SHEET

Cannula size	25G
Cannula guard	Yes
Classification	Dispersive
Composition	30 mg/ml (3% sodium hyaluronate)
Volume	0.85 ml
Viscosity	75,000 mPas
Molecular weight	1,000,000 Daltons
Storage temperature	+2°C to +30°C
Substance	Fermented HA; animal-free
SKU	10301012

- Refrigeration not required¹
- Animal-free Hyaluronic Acid¹
- 98% of surgeons felt satisfied with adherence of Healon EndoCoat™ PRO OVD to the corneal endothelium^{2*}
- 99% of surgeons felt satisfied with the transparency and visibility of Healon EndoCoat™ PRO OVD to the corneal endothelium^{2*}

*Based on responses of 11 surgeons who answered "Completely satisfied" or "Mostly satisfied" following use in multi-center Clinical Trial (n=131)

REFERENCES:

1. HEALON EndoCoat™ PRO IFU, 52-0629-00E, current revision
2. 2025REF5178 Clinical Trial Questionnaire Results

The promoted product is a medical device. For healthcare professionals only. Please reference the Instructions for Use for a complete list of Indications and Important Safety Information and contact our specialists in case of any question.

© Johnson & Johnson and its affiliates 2026.

2025PP16703

Johnson & Johnson

 **Healon™**
FAMILY of OVDs